


Assessing the Cardiotoxicity of 31 Compounds Using a Multiplexed Kinetic Image Cytometry® (KIC®)-based Assay: Harnessing the Predictive Power of Human iPSC-Cardiomyocytes for Early Drug Development



Ranor C. B. Basa^{*†a}, Randall S. Ingermanson^a, Filiberto Catalan-Perez^a, Ricardo Serrano^b, Ariel Wang^a,

Alyson Smith^a, Jeffrey M. Hilton^a, Patrick M. McDonough^a, Cherie Handley^a, Lilian Harrison^a, Ameena Pascua^a,

Mark Mercola^b, Jeffrey H. Price^{†a,c}

^a Vala Sciences, Inc., 6370 Nancy Ridge Dr., Suite 106, San Diego, CA, USA


^b Stanford Cardiovascular Institute, Stanford University, Stanford, CA, USA

^c Scintillon Institute, San Diego, CA, USA


*📧 [†]✉️ ^arbasa@valasciences.com, jprice@valasciences.com

Automated KIC® Measurements: Examples from Assay Development


Single-Cell CaT/AP




Frequency




Peak Duration




EADs/DADs



Contractility



Mito Membrane Potential




• Traditional live-cell fluorescence microscopy (single image, non-kinetic analysis)

Statistical Analysis


- **Bar plots:** bars and error lines represent mean±S.D. (this section only)
- **Significance testing:** Ordinary one-way ANOVA followed by Dunnett's post-hoc multiple comparisons test (*p<0.05, **p<0.01, ***p<0.001) for N=3 (generally) using GraphPad Prism 9

CyteSeer® automatically computes >200 measurements per optical channel!


Methods: KIC® Cardiotoxicity Assay



Methods: CyteSeer® Analysis



Waveform classification and predicted ECG correlations from functional dose-response relationships




Results Summary: 31-Compound Pilot Study

Drug #	Drug Name	Clinical Cardiac Effect*	Mechanism of Action/Cardiac Target Class	Cardiac Target(s)	Target (AG)onist or (ANT)agonist	Measurement(s) of Cardiac Effect† [Dose (μM)]				Predicted Human Outcome?
						AP	SR	CON	MT	
1	Milrinone	HF- (HF+, Arr+)	SR	PDE3	ANT	[50] ▲APD75, ▲Ampl, ▲BPM, ▲UVel, ▲DVel	[50] ▲DVel	[50] ▲DVel	ns	Yes
2	Flecainide	Arr-	SR	Nav1.5 (I _{Na}), RyR2	ANT, ANT	[9.49] Q	[30] ▼Ampl, ▼UVel, ▼DVel	[0.095] ▼MTD75; [9.49] ▼Ampl, Q	ns	Yes
3	Istaroxime	HF-	SR	SERCA2, Na ⁺ /K ⁺ ATPase	AG, ANT	[0.0125] ▼APD75; [0.025] ▲Tri	[0.0125] ▼CTD75	[0.025] ▼MTD75; [0.05] ▲UVel, ▲DVel; [0.1] ▲Ampl	ns	Yes
4	Cisapride	Arr+	AP	hERG (I _{Kr})	ANT	[0.125] ▲APD75	[0.125] ▼DVel	[0.25] ▲MTD75	ns	Yes
5	Diltiazem	Arr-	AP	Cav1.2 (I _{CaL} /I _{CaT})	ANT	[0.625] ▼BPM, ▼APD75	[1.25] ▼Ampl, ▼CTD75	[1.25] ▼MTD75	[2.5] ▼MtMP	Yes
6	Propranolol (racemic mixture)	Arr-	AP	β ₁ -AR, R(+)-Nav1.5 (I _{Na})	ANT, ANT	[3.125] ▼BPM	[3.125] ▼UVel, ▼DVel	[12.5] ▼DVel	[25] ▼MtMP	Yes
7	Mexiletine	Arr-	AP	Nav1.5 (I _{Na}), hERG (I _{Kr})	ANT	[0.32] ▲UVel	[3.125] ▼DVel	[6.25] ▲MTD75	[100] ▼MtMP	Yes
8	Sotalol	Arr+	AP	hERG (I _{Kr})	ANT	[3.25] ▼BPM	[6.25] ▼DVel	[3.16] ▼MTD75	ns	Yes
9	Ivabradine	HF-	AP	HCN (I _h), hERG (I _{Kr})	ANT, ANT	[1.25] ▲Tri	[0.625] ▼DVel	[1.25] ▼Peak#	[10] ▼MtMP	Yes
10	Ibutilide	Arr- / Arr+	AP	hERG (I _{Kr})	ANT	[0.0156] ▲APD75	[0.0156] ▲CTD75	[0.0156] ▼Peak#	ns	Yes
11	Quinidine	Arr-	AP	hERG (I _{Kr})	ANT	[0.156] ▼DVel	[0.312] ▼DVel, ▼BPM	[0.312] ▼DVel	ns	Yes
12	Dofetilide	Arr+ / Arr-	AP	hERG (I _{Kr})	ANT	[0.0156] ▲Tri	[0.0156] ▲CTD75	[0.0156] ▲MTD75	[0.5] ▼MtMP	Yes
13	Isoproterenol	Arr-	AP	β-AR, K _v 7.1 (I _{Kr})	AG, AG	[0.01] ▼Tri, ▲DVel	[0.01] ▲Ampl, ▲DVel	[0.03] ▲DVel	[10] ▼MtMP	Yes
14	Trimetazidine	HF-	CON / MT	Mt LC 3-KAT	ANT	[50] ▲Tri	[50] ▼DVel	ns	ns	No
15	Levosimendan	HF-	CON	Ca ²⁺ -TnC	AG	[50] ▲APD75	[25] ▲Ampl, ▲CTD75	[100] ▲Ampl	[100] ▲MtMP	Yes
16	Omecamtiv	HF- (failed trial)	CON	Myosin	AG	[2.5] ▼Tri, ▲APD75	ns	[0.625] ▼UVel, ▼DVel; [10] Q	ns	Yes
17	Mavacamten	HF-	CON	Myosin	ANT	[0.095] ▲APD75	[5] ▼CTD75	[0.095] ▼Peak#	ns	Yes
18	Vandetanib	Arr+, HF+	All	hERG (I _{Kr})	ANT	[1.25] ▲APD75	[1.25] ▼BPM	[0.625] ▲MTD75	ns	Yes
19	Sorafenib	HF+	All	SERCA2 / cytotoxic	ANT	[0.00032] ▲UVel	[0.0316] ▼Ampl	[0.00316] ▼Ampl	ns	Yes
20	Regorafenib	HF+	All	I _{NaL}	AG	[31.6] ▼BPM	[31.6] ▼UVel	[50] ▼Peak#	ns	Yes
21	Carboplatin	HF+	All	I _{CaL}	AG	[12.5] ▲APD75	[12.5] ▲CTD75	[6.25] ▲MTD75	ns	Yes
22	Ranolazine	Arr-	SR / AP	RyR2 / I _{Na}	ANT	[25] ▲APD75, ▲Tri	[6.25] UVel	[25] ▲MTD75	ns	Yes
23	Veratridine	Arr+ / Q?	AP	Nav1.5, Nav1.8 (I _{Na})	ANT	[0.312] ▲APD75	[0.312] ▲CTD75	[0.156] ▼Peak#	ns	Yes
24	Thapsigargin	Q	SR	SERCA	ANT	[0.625] ▼APD75	[0.312] ▼Ampl	[0.625] ▼Ampl, ▼MTD75	ns	Yes
25	Caffeine	Neg	Neg	-	-	ns	ns	ns	ns	Yes
26	Loratadine	Neg	Neg	-	-	[0.2] ▲APD75	[0.2] ▲CTD75	ns	ns	Yes
27	Acetaminophen	Neg	Neg	-	-	ns	ns	ns	ns	Yes
28	Doxycycline	Neg	Neg	-	-	ns	ns	ns	ns	Yes
29	Ciprofloxacin	Neg	Neg	-	-	ns	ns	ns	ns	Yes
30	Oseltamivir	Neg	Neg	-	-	ns	ns	ns	ns	Yes
31	Aspirin	Neg	Neg	-	-	ns	ns	ns	ns	Yes

* **Clinical Cardiac Effect:** Arr=arrhythmia (+, causes; -, treats); HF=heart failure (+, causes; -, treats); Q=arrest (quiescence); Neg=negative control/no clinical cardiac effect
† **Measurement(s) of Cardiac Effect:** [X]=lowest dose at which significant changes were observed, in μM; ▲=increase/▼=decrease; AP=action potential/membrane voltage; SR=sarcoplasmic reticulum/Ca²⁺ transient; CON=contraction (APD/CTD/MTD75=Action Potential/Ca²⁺ Transient/Motion [contractile] Transient Duration, 75% down from peak; UVel/DVel=max rate of peak upstroke/downstroke; Ampl=peak amplitude; Tri=triangulation; BPM=beats/min; Peak#=number of contractile peaks); MT=mitochondria (MtMP = mitochondrial membrane potential); ns=no significant changes observed

prophECG: “Predicting” the Safety of Vandetanib



- o KIC® of hiPSC-CMs provides **sensitive, data-rich readouts of multiplexed voltage, calcium, contraction, and mitochondrial membrane potential.**
- o *In vitro* KIC® measurements **correspond to *in vivo* clinical readouts (e.g., ECG);** large KIC® data sets can be used to **predict compound effects in humans, thus accelerating early-stage drug discovery.**
- o In the next year, we **will screen 300+ drugs** with known cardiotoxicity risk profiles (e.g., from databases like DICTrank^{13,14}), as well as re-analyze existing KIC® data (e.g., from Pfeiffer *et al.*⁷) with these newer methods.
- o We continue to develop and refine computational strategies (including **classical AI/ML**) to provide further insight into the complex interplay among all the different KIC® measurements in these data for **improved, explainable predictivity and more clinically-relevant (and potentially novel) classifications.**

References

1. Clark, M., & T. Steger-Hartmann. A big data approach to the concordance of the toxicity of pharmaceuticals in animals and humans. *Regul Toxicol Pharmacol* 96, 94-105 (2018). PMID: 29739468.

2. Valentin, J.P., et al. The Challenges of Predicting Drug-Induced QT Prolongation in Humans. *Toxicol Sci* 147, 3-24 (2022). PMID: 35148601.

3. Maron, L.M., et al. The automated gains pg an effective early cardiovascular screening and lead optimization model. *J Pharmacol Toxicol Methods* 68, 137-149 (2013). PMID: 23600010.

4. Dewart, L., et al. The concordance between nonclinical and phase I clinical cardiovascular safety data from a cross-company data sharing initiative. *Toxicol Sci* 142, 47-59 (2014). PMID: 24546661.

5. Bui, S., et al. Predicting the Clinical Translation of Pharmacokinetic Effects on Cardiovascular Safety Pharmacology Studies. *Toxicol Sci* 146, 272-281 (2019). PMID: 30736991.

6. Lu, H.K., et al. High-Throughput Measurement of Ca²⁺ Dynamics in Human Stem Cell-Derived Cardiomyocytes by Fluorescence Microscopy. *J Cardiovasc Pharmacol* 54, 205-210 (2012). PMID: 22630003.

7. Pfeiffer, E.R., et al. Specific prediction of clinical QT prolongation by kinetic image cytometry in human stem cell-derived cardiomyocytes. *J Pharmacol Toxicol Methods* 87, 253-271 (2016). PMID: 26908051.

8. Pfeiffer, E.R., et al. High-throughput screening of cytotoxicity using kinetic image cytometry in human stem cell-derived cardiomyocytes. *J Pharmacol Toxicol Methods* 87, 253-271 (2016). PMID: 26908051.

9. Conroy, F., et al. High-throughput measurement of Ca²⁺ dynamics for drug risk assessment in human stem cell-derived cardiomyocytes by kinetic image cytometry. *J Pharmacol Toxicol Methods* 86, 246-256 (2015). PMID: 25925523.

10. Sefton, L., et al. A deep learning platform to assess drug proarrhythmia risk. *Cell Syst* 30, 96-99 (2020). PMID: 32605085.

11. Pfeiffer, E.R., et al. Predicting the clinical translation of pharmacokinetic effects on cardiovascular safety pharmacology studies. *Toxicol Sci* 146, 272-281 (2019). PMID: 30736991.

12. Pfeiffer, E.R., et al. Predicting the clinical translation of pharmacokinetic effects on cardiovascular safety pharmacology studies. *Toxicol Sci* 146, 272-281 (2019). PMID: 30736991.

13. Qiu, Y., et al. The largest reference list of 1218 human drug-induced cardiotoxicity using FDA labeling. *Drug Safety Toxicol* 11, 1037-1051 (2023). PMID: 37714400.

14. Sefton, L., et al. Insights into Drug-Induced Cardiotoxicity from Biological and Chemical Data: The Public Classification for FDA Drug-Induced Cardiotoxicity Risk. *J Chem Inf Model* 64, 4, 1175-1186 (2024). PMID: 38308051.

15. We would like to thank the following individuals for their technical assistance and/or advice: Ramona Aguinaga (Vala Sciences, Inc.), Dries A. W. Foppen and Michelle Xu (Stanford University).

16. This work was supported by a grant from NIH/NHLBI (1R01HL158510, PI: Jeffrey H. Price).

Acknowledgments

Now seeking pharma partners for validating our approach!